




RMK700 Keratometer


Exquisite workmanship, high quality design, smooth movement together with pure electric lifting, patient will experience the relaxing measurement process. For application of Hartmann technology, every data is the result of calculation and comparison overing 20 times. Both accuracy and rapid of measurement can meet all requirement of the patient.




A wide variety of measurement mode & function




REF Mode




KER Mode



K&R Mode




SIZE Mode




CLBC Mode



Concise menu interface



MAP Mode



One step locking switch of body head




Convenient and sample lifting forehead rest



Electric lifting chin rest



Perfect optical head and foggy system are much more accurate



Clear foggy interface make the measuring possess more comfortable



Integrated printer is more convenient and faster

RMK700

With function of REF, KER, K&R, CLBC, SIZE, RMK-700 can meet all of your measurement requirement.

Convenient trying on mode, patient can experience the measurement result on the machine. Omitting all the trouble of trying again and again.

Z-map measurement mode is much more vivid than traditional cornea measuring. By the result of outputting to PC, the patient can understand all the details of his vision.



Performance Parameter Refractometer

Vertex Distance (VD): 0.0 mm, 12.00 mm, 13.5 mm, 15.00mm

Sphere Power (SPH): -20.00m-1 ~ +20.00m-1 (VD=12mm, 0.12 m-1/0.25 m-1)

Cylinder Power (CYL): 0.00 m-1~+10.00 m-1 (0.12m-1/0.25m-1)

Axis (AX): 1°~ 180° (1°)

Cylinder Form: -, +, MIX

Minimum Pupil Diameter: 2.0mm

Keratometer (only apply to model with KER function)

Radius of Curvature: 6.5 mm ~ 9.4 mm (Increment: 0.02mm)

Others

Size: 506mm×269mm×477mm

Monitor: LCD Color Monitor

Supore

RMK700 Keratometer



Supore®

Company: Shanghai Supore Instruments Co., Ltd.
Add: No.800 Yeji RD Shanghai China

Tel: +86-21-36160506
Fax: +86-21-36162800-888

E-mail: sales@supore.com
Http:www.supore.com