

ACP-60S

AUTO CHART PROJECTOR



Supore

- Providing several visual target and visual function icons
- Small volume
- Easy to operate



ACP-60S

www.supore.com

Liquid Crystal Eye Chart

- More clear images
- Switching visual targets smoothly and quickly
- Space-saving wall-hanging design
- Thinner and more aesthetic design

The input power supply parameters:
AC100-240V/50HZ/60HZ
Standard power consumption: About 45W
Power consumption: Less than 2W
Maximum screen resolution: 1440 x 900
Screen dot spacing: 0.248mm
Screen type: TFT

Screen ratio: 16:9 or 16:10 (slight difference is acceptable according to the actual batch)
Screen size: 21.5
Maximum brightness: 300cd/m
Maximum contrast: 2000:1
Signal response time: 5ms
Boundary dimension: 515mm*320mm*50mm

Use Distance: 1-7 m
Visual acuity record method:
The decimal English scores recording
Remote-controlled distance: No less than 8 m
Installation method: Bracket can be installed and hung on the wall