





LIDYA

OPHTHALMIC EXAMINATION UNIT



MABEL
mabelas.com.tr



LIDYA

OPHTHALMIC EXAMINATION UNIT

 **MABEL**
mabelas.com.tr





	BASIC TECHNICAL CAPACITY	STANDARD	OPTIONAL
1	BASE PLATE, MAIN BODY, SLIDING TABLE (FIXED HEIGHT) FOR 2 INSTRUMENTS	*	
2	MOTORIZED PATIENT CHAIR (UP & DOWN/24 CM.) BACK IS FIXED	*	
3	COLUMN FOR CHART/PHR. ARM. & READ LAMP	*	
4	READ (HEAD) LAMP	*	
5	DRAWER FOR TRIAL LENS SET		*
6	MECHANICAL BREAK	*	
7	TIPUP ARM REST	*	
8	UNIVERSAL CHINREST PLATE FOR SLITLAMPS	*	
9	OUTPUTS FOR SLITLAMPS LED & HALLOGEN (5V/6V/12V/18V)	*	
10	SLIDING TABLE 90 DEGREE ROTATION & PATIENT CHAIR 360 DEGREE FUNCTIONAL MOVEMENT WITH LOCK	*	
11	RIGHT & LEFT HAND SELECTION	*	
12	SECURITY SWITCH UNDER SLIDING TABLE TO PROTECT PATIENT	*	
13	220 V OUTPUT FOR EXTRA DEVICES & SLEEPING MODE	*	
14	SLIDING TABLE MOTORIZED UP & DOWN MOVEMENT (14 CM.) WITH FOOTSWITCH		*
15	RECLINER PATIENT CHAIR		*
16	CHART PROJECTOR HOLDER		*
17	ELECTROMAGNETIC BREAK (TO LOCK THE TABLE ALL POSITIONS) "WITH MECHANICAL BRAKE TOGETHER"		*
18	PHOROPTER ARMS (MOTORIZED OR MANUEL TYPES)		*
19	TIP UP FOOTREST		*
20	DESK FOR PHOR. CONTR. PANEL OR LENSMETER		*
21	WHEELCHAIR ACCESSABLE BASE PLATE		*
22	IT CAN BE ASSEMBLED FOR BOTH RIGHT AND LEFT-HAND USE (STANDARD)	*	*

LIDYA OPHTHALMIC EXAMINATION UNIT TECHNICAL DATA

	Standard	Min. & Max.	Width & Length	OUTPUT	210-240VOLT AC
Patient chair height		520 mm 760 mm		Read Lamp	12-24VOLT/10WATT
Sliding table(Fixed) height	825 mm			Fuses	5A
Sliding Table(Up&Down Motorized "optional") height		825 mm 965 mm			
Sliding table dimensions			425 mm 940 mm		
Drawer for TLS dimensions			468 mm 680 mm		
Patient Chair Load		Max 180 kg			

info@mabelas.com.tr

+90 312 256 33 44


Ostim OSB Mah. 1220. Cad.
No:35 Yenimahalle/ANKARA

Social Media



mabelmakina