

Pi ADVANCE

OPERATION TABLE



ADVANCE

PI ADVANCE

OPERATION TABLE

Optimal Comfort,
Maximum Precision.



Pi ADVANCE

OPERATION TABLE



Color Options:



Jet Black Ruby Wine Forest

PI ADVANCE

OPERATION TABLE



PI ADVANCED MODEL EYE OPERATION TABLE TECHNICAL DATA

	Standard	Minimum	Maximum
Table Height		540 mm	850 mm
Headrest (diameter)	200 mm		
Backrest	610 mm x 640 mm (h)		
Sitting Place	570 mm x 640 mm (h)		
Legrest	460 mm x 610 mm (h)		
Length		1680 mm	2080 mm
Max. Width	720 mm		
Max. Lifting Capacity	4000 N (400 kg)		
Safe Lifting Capacity	300 kg		
Wheel (diameter)	100 mm		
Backrest Angle Degree	180/115		
Legrest Angle	180/110		

TECHNICAL SPECIFICATIONS
BATTERY TYPE: LI-ION, LITHIUM
BATTERY CAPACITY: 2.5/4.5/10Ah (LITHIUM)
VOLTAGE: 24V DC, 2.9/5A
CHARGE VOLTAGE: 29~45V DC
IP RATING: IP66
CERTIFICATE: IEC60601-1, ES606061-1
OPERATIONAL TEMPERATURE RANGE: +5°C~+45°C
CHARGE INDICATION (FOR CHARGE/DISCHARGE STATUS)
BUILT-IN BUZZER
STANDBY POWER CONSUMPTION
DETACHABLE BATTERY PACK, BATTERY CAPACITY INDICATOR

madde 16

madde 17

madde 15

madde 19






madde 18

PI ADVANCE

OPERATION TABLE

MABEL
mabelas.com.tr

Technical Dimensions:



MABEL
mabelas.com.tr

PI ESSENTIALS

OPERATION TABLE



ESSENTIALS

PI ESSENTIALS

OPERATION TABLE

Designed for Surgeons,
Perfected for Patients.



Color Options:




Jet Black Ruby Wine Forest Green

PI ESSENTIALS

OPERATION TABLE

MABEL
mabelas.com.tr

Technical Dimensions:



PI ESSENTIALS MODEL EYE OPERATION TABLE TECHNICAL DATA


	Standard	Minimum	Maximum
Table Height		540 mm	850 mm
Headrest (diameter)	200 mm		
Backrest	865 mm x 640 mm (h)		
Legrest	875 mm x 640 mm (h)		
Max. Width	720 mm		
Max. Lifting Capacity	4000 N (400 kg)		
Safe Lifting Capacity	300 kg		
Wheel (diameter)	100 mm		

BATTERY TECHNICAL SPECIFICATIONS

BATTERY TYPE: LEAD ACID, LITHIUM
CAPACITY: 2.9/5Ah (LEAD ACID); 2.5/4.5/10Ah (LITHIUM)
DISCHARGING VOLTAGE: 24V DC, 2.9/5A
CHARGING VOLTAGE: 29~45V DC
COLOR: GREY
IP RATING: IP66
CERTIFICATE: IEC60601-1, ES606061-1
OPERATIONAL TEMPERATURE RANGE: +5°C~+45°C
LED INDICATION (FOR CHARGE/DISCHARGE STATUS)
WITH BUILT-IN BUZZER
LOW STANDBY POWER CONSUMPTION
DETACHABLE BATTERY PACK, BATTERY CAPACITY INDICATOR

info@mabelas.com.tr

+90 312 256 33 44



Ostım OSB Mah. 1220. Cad.
No:35 Yenimahalle/ANKARA

Social
Media

Instagram Facebook LinkedIn mabelmakina