

TETA

SURGEON CHAIR



TETA

SURGEON CHAIR

MABEL
mabelas.com.tr



TETA

SURGEON CHAIR

 **MABEL**
mabelas.com.tr



TECHNICAL DATA

		STANDARD	MINIMUM	MAXIMUM
1	BACKREST	380 x 390 mm (h)		
2	SITTING PLACE	490 x 450 mm (h)		
3	LENGTH		440 mm	680 mm
4	MAX. WIDTH	720 mm		
5	MAX. LIFTING CAPACITY	2000 N (200 kg)		
6	SAFE LIFTING CAPACITY	180 kg		
7	WHEEL (DIAMETER)	100 mm		
8	BACREST ANGLE DEGREE	85 / 115		



info@mabelas.com.tr

+90 312 256 33 44



Ostim OSB Mah. 1220. Cad.
No:35 Yenimahalle/ANKARA

Social
Media

Instagram Facebook LinkedIn mabelmakina